

October 2012

Nerve Stimulator with High Reliability for Assistance in Regional Anesthesia

By **Carlos Alexandre Ferri**, Master's candidate at University of Campinas in Brazil with special acknowledgement to **Professor Antonio A. F. Quevedo**, Assistant Professor at the Department of Biomedical Engineering, Faculty of Electrical and Computer Engineering



freescale.com/medical

1 Introduction

Plexus blockade, or regional anesthesia, is the technique used in minor surgery to eliminate pain, such as in surgery of the hand, arm, elbows, feet, knees, legs and so on. Its preference over general anesthesia is due to a shorter stay in recovery room, minimized systemic adverse effects, decreased incidence of nausea and vomiting and eliminated need for intubation (Cangiani et al, 2006).

In the last decades, the use of nerve stimulators for identifying plexus in assisting regional anesthesia has been a technique of great importance for the identification of nerves, consolidating highly beneficial to the patient, because (Bollini et al, 2006, Regatieri, 2002):

- Reduces risk of injury to the nerve;
- Allows a precise identification of the nerve plexus;
- Possible reduction in the dose of anesthetic blockade for providing best quality;
- Reduction of collateral effects from excess anesthetic.

These benefits are possible due to the ability of the nerve stimulator to estimate the distance between the needle tip and the nerve. Therefore, one can identify the best place to infiltrate the local anesthetic solution (Dalrymple et al, 2006).

1.2 Use of the Nerve Stimulator to Aid Plexus Location

Nerve stimulation occurs only when current intensity is applied by the needle electrode and the tip is sufficiently close the nerve. Consequently, a muscle contraction becomes evident (Andrés et al, 2001).

The procedure for locating plexus and the correct placement of the needle by use of the nerve stimulator begins with a standard routine for regional anesthesia, followed by adjustment of the initial current stimulator for 1 to 2 mA and the duration of the stimulus for 100us, with a frequency of 2 Hz. The needle is inserted until a contraction or muscle spasm is identified visually. After identification

of the contraction, current amplitude is decreased until the contraction disappears. Introduction of the needle continues until a new contraction appears. This procedure is repeated until the applied current is 0.3 mA and a muscle contraction is visible. At this point, the needle tip is from 1 to 2 mm distance from motor nerve and the injection of local anesthetic solution should provide a satisfactory blockade (Sardesai et al, 2009).

1.2 Features Standard Equipment

Much of the equipment that is available on the market designed to aid the regional anesthetic procedures does not have the minimum requirements recommended for a good neurostimulator (Hadzic et al, 2003; Jochum et al, 2006). This shortcoming in the equipment can increase the risk to the patient, since the minimum requirements ensure its use in the identification process of the nerve plexus with the precision necessary to achieve a satisfactory blockade.

Constant current: The ability of the neurostimulator to provide constant current is of critical importance for achieving plexus blockade. Tissue impedance may vary between different individuals. These variations may be caused by diseases like diabetes or renal problems, and also higher impedance obstacles such as bones. Thus, the equipment must be able to provide a constant current to a load variation between 1kΩ and 10kΩ (Barthram, 1997; Jochum et al, 2006).

Duration of pulses: The duration of the stimuli should correspond to chronaxy of motor neurons to ensure that Aα sensory neurons are not stimulated and hence muscle contraction occurs in the absence of pain. Thus, the

stimulator should have the following options for stimulus duration: 50μs, 100μs, and 200μs (Barthram, 1997; Jochum et al, 2006).

Accuracy: Accurate control of stimulation parameters is mandatory for the correct placement of the needle. This means that the current delivered by the stimulator must be equal to the expected current, the morphology of the pulses must be rectangular, the pulse width variation should be less than five percent of the expected value, and rise and fall times should be less than 3μs (Hadzic et al, 2003; Jochum et al, 2006).


Frequency Stimulation: The frequency of stimulation allows control of the speed at which the nerve is located. The frequency should be selected to 1 Hz or 2 Hz, the most usual frequency of 2 Hz (Barthram, 1997; Jochum et al, 2006).

Functionality: The stimulator must provide maximum safety for the patient and provide the user with pertinent information about the operation of the system, such as visual and audible warning at the time of stimulation; battery information; high impedance warning; and indication of selected parameters (Regatieri, 2002; Jochum et al, 2006).

2 Method

The stimulator developed in this project can be summarized by the block diagram shown in Figure 1. The system is organized into modules to facilitate identification of possible sources of noise and failures. As central control processes, we opted to use the Freescale development board DEMO9S08LC60. This board is responsible for the communication interface between the stimulator and the PC (that is isolated from

Figure 1: Block Diagram of Operation of the Stimulator



the power grid, powered by battery), control of stimulation parameters and graphical user interface via an LCD display. Modules 4 and 5 will be added in the future as a feedback mechanism to the neurostimulator, allowing the current amplitude settings to be performed automatically. This control will give more autonomy to the anesthesiologist.

2.1 Module 1: Supply

This equipment has to be portable for convenience. Therefore, its power supply is derived from two 9V batteries in series, providing 18V. Whereas the maximum load impedance does not exceed 10kΩ, and stimulus current is never beyond 5 mA, the voltage compensation will be a maximum of 50V. We decided to use a charge pump converter, Dickson type.

The output voltage of this circuit can be calculated by the following formula:

(1)

$$V_{out} = (N + 1)V_{DD} - N \frac{I_L}{C}$$

Where N is the number of stages, VDD is the supply voltage, IL is the current output to the load, C capacitance and Vout is the output voltage.

Considering the need for a voltage of 50V, and that the voltage supplied by the battery may vary between 18V and 15V during its useful life, the converter should have three


stages (n = 3). Thus, there is a Vout ranging between 60V and 72V. A voltage regulator coupled to the output of the circuit ensures that the output voltage is stabilized at 50V. This regulator consists of a bipolar transistor biased by a zener diode (Vz = 52V) connected at its base. A negative voltage of -5V was obtained through an inverter to a charge pump stage with its output set to -5V by a bipolar regulator. The voltages of 5V and 12V were determined using bipolar regulators.

2.2 Module 2: Pulse Generator and Digital Controller

This module is aimed at controlling the duration of the stimulus and amplitude current digitally. An open-collector comparator was used to control the width of the pulses (CLK), which is enabled by the MCU. To control current amplitude, a DAC (digital-to-analog converter), connected to a pull-up resistor (R) which is connected to the output of the comparator controls the voltage level (VDAC) necessary for generating the current required. A constant current source is accomplished by a current voltage converter comprising the amplifier IC2, transistor Q1 and RREF (Figure 2). Thus, the current generated is given by:

$$I_{REF} = \frac{V_{DAC}}{R_{REF}} \times CLK$$

Figure 2: Schematic of Control and Generator Stimulus



2.3 Module 3: Current Mirror

As human body impedance may vary due to several factors, a system that provides a constant current independently of load change was implemented. The characteristic of a current mirror fits this application, allowing a constant current, invariant with the load. Thus, the current mirror must have high output impedance so that the current does not undergo changes with the load. Wilson type was chosen because it has a high output impedance and large output voltage excursion.


The current measurements to check the linearity of the circuit were made with a multimeter connected in series with output impedance.

2.4 Development Board

The DEMO9S08LC60 development system was used for digital control of the other modules. This board integrates the Freescale MCU MC9S08LC60, plus a built-in LCD display, RS232 interface and buzzer. This board provides all the tools needed to develop this nerve stimulator because it is able to control all stimulation parameters, such as frequency and duration of the stimuli by triggering timers, current amplitude by updating the DAC through SPI (serial peripheral interface), communication with PC through RS232 interface, plus visual and audible warnings through LCD and buzzer, respectively. The flow diagram shown in Figure 3 illustrates the sequence of work performed by the MCU.

Subsequently, the equipment will receive a module for capturing electromyographic signals (EMG), whose purpose is to provide a feedback for process automation of current level control. Thus, the AD converter of the MCU will be used. Complementary to this process, the use of a Freescale accelerometer (MMA84530) attached to a selected muscle enables data collection from a muscle contraction evoked by stimulation. These data enhance the feedback process in the automatic control of current levels.

Figure 3: Sequence of Work Performed by the MCU



2.5 PC Interface

All settings to the nerve stimulation and location are controlled through a graphical user interface developed for Windows® environment in Microsoft Visual Basic®. This interface allows the user to select stimulus duration, frequency and current amplitude with ASCII commands sent to the MCU through RS232 serial interface at a transmission rate of 115 kbps.

3 Results


A very important factor for electrical stimulation with invasive electrodes is low-noise generation from electronics. The major source of noise in this system originates in voltage converters. Thus, we observed the noise generated at the outputs of the converters (Table 1).

Table 1: Noise Level in the System

	+5V	-5V	+12V	+50V
Noise (mV)	-	-	-	0,1

The output voltage of the Dickson topology remained constant until the output current was 4 mA, during which it suffered a significant decrease. This decline was expected, as seen in Equation 1. However, this variation in output voltage is not an issue to the system since the current pulses are of short duration (less than 200 μs), which causes a voltage drop lower than 200 mV when applying a current pulse of 5 mA with duration of 200 μs in a load of 8 kΩ. It also exhibited rapid recovery.

Figure 4: Linearity For the Wilson Model



The pulses of current generated by the system presented rise time of less than 0.4 μs , 0.25 μs for falling time; width of the pulses with variation of less than one percent.

The current mirror presented excellent results in maintaining a constant current, independent of the load, as seen in Figure 4. This means that the circuit is able to generate and maintain the actual current equal to the expected current within the load range. The output impedance of the current mirror was measured indirectly by means of voltage variation observed for two different loads when the same current is selected. Thus, the calculated output impedance is 1.7 $M\Omega$.

4 Discussion

Electrical noise exceeding 10 mV may trigger undesired stimulation. The boost converter was chosen first for this system. However, this circuit showed high noise level, in the order of 100 mV, precluding its use. The charge pump converter and Dickson topology presented

satisfactory results, although its efficiency is not high (75 percent), when compared to other conversion techniques, which have efficiency greater than 80 percent.

To make sure the reference current was linear with high accuracy, the DAC should provide excellent linearity. The converter used for this project was adequate, with a coefficient of linearity (R^2) equal to 1. This provides assurance that the actual current is equal to the expected current.

The stimuli generator responded appropriately, with rise and fall times well below the maximum expected. This ensures that the stimulus format is closest to the rectangular, as desired.

The graph in Figure 4 reflects that the system has excellent linearity. Thus, the load current (i.e., the current in the patient) will be within the range of accuracy to the stimulation purpose of locating plexus.

5 Conclusion

The development equipment was of excellent quality, and all of the features needed for proper equipment were met. The low variation in stimulation parameters such as pulse width and morphology was analyzed with excellent correlation between the expected and the real current. These are critical factors for selective stimulation of motor fibers $A\alpha$. During induction, the patient submitted to the use of the equipment to avoid discomfort and pain. The accuracy in the parameters is what ensures an excellent plexus location.

The equipment works with a PC interface (in this case, a battery-powered notebook computer that was isolated from the power grid), which allows control of frequency, width and amplitude of the stimulation current. Moreover, this communication will allow storing the information obtained through EMG and 3-axial acceleration of muscle contraction to be evoked by electrical stimulation. This information will be used for automating the neurostimulator, and also for assessing the patient's motor response to local anesthetics, as well as a pre-test of the quality of regional blockade.

After field testing, we intend to change the MCU so the final equipment will be more portable, with less power consumption. Our probable choice will be the Kinetis K53 MCU which has these embedded features: LCD controller, digital-to-analog converter, touch sensing interface and an op-amp that can be used as a PGA. With all of these features, the final prototype will demand less external components.

6 References

André's, J. de; Sala-Blanch, X. "Peripheral Nerve Stimulation in Practice of Brachial Plexus Anesthesia: A review." *Regional Anesthesia and Pain Medicine*, Vol 26, no 5, pg 478-483, 2001.

Barthram, C. N. "Nerve stimulators for nerve location - are they all the same.," *Anaesthesia*, vol. 52, pg. 761-764, 1997.

Bollini, C. A.; Cacheiro, . "Peripheral nerve stimulation". *Techniques in Regional Anesthesia and Pain Management*, Vol 10, No 3, July 2006.

Cangiani, L. M.; Posso, I. P.; Braga, G. M.; Nogue, C. S. *Tratado de Anestesiologia. SAESP/ Organização. Sociedade de Anestesiologia do Estado de São Paulo*. 6. ed. 2006.

Dalrymple, P.; Chelliah, S. "Electrical nerve locators". *Continuing Education in Anaesthesia, Critical Care & Pain*. Vol 6, No 1, pg 32-36, 2006.

Hadzic, A.; Volka, J.; Hadzic, N.; Thys, D. M.; Santos, A. C. "Nerve Stimulators Used for Peripheral nerve blocks vary in their electrical characteristics," *Anesthesiology*, no. 98, pg. 969-74, 2003.

Jochum, D.; Iohom, G.; Diarra, D. P.; Loughnane, F.; Dupré, L. J; Bouaziz, H. "An objective assessment of nerve stimulators used for peripheral nerve blockade," *Anaesthesia*, vol. 61, pg. 557-564, 2006.

Palumbo, G; Pappalardo, D. "Charge Pump Circuits: An Overview and Design Strategies and Topologies". *IEEE Circuits and Systems Magazine*. First Quarter, 2010.

Regatieri, F. L. F. "Uso do estimulador de nervos em bloqueios periféricos." *Revista da Sociedade de Anestesiologia do Estado de São Paulo*, vol. 4, pg. 18-20, 2002.

Sardesai, A. M.; Iyer, Uma. "Nerve stimulation for Peripheral nerve Blockade." *Anaesthesia Tutorial of the Week 149*. Disponível em: <http://totw.anaesthesiologists.org>, 2009.

How to Reach Us:

Home Page:

freescale.com

e-mail:

support@freescale.com

USA/Europe or Locations Not Listed:

Freescale Semiconductor
Technical Information Center, CH370
1300 N. Alma School Road
Chandler, Arizona 85224
1-800-521-6274
480-768-2130
support@freescale.com

Europe, Middle East, and Africa:

Freescale Halbleiter Deutschland GmbH
Technical Information Center
Schatzbogen 7
81829 Muenchen, Germany
+44 1296 380 456 (English)
+46 8 52200080 (English)
+49 89 92103 559 (German)
+33 1 69 35 48 48 (French)
support@freescale.com

Japan:

Freescale Semiconductor Japan Ltd.
Headquarters
ARCO Tower 15F
1-8-1, Shimo-Meguro, Meguro-ku,
Tokyo 153-0064, Japan
0120 191014
+81 3 5437 9125
support.japan@freescale.com

Asia/Pacific:

Freescale Semiconductor Hong Kong Ltd.
Technical Information Center
2 Dai King Street
Tai Po Industrial Estate,
Tai Po, N.T., Hong Kong
+800 2666 8080
support.asia@freescale.com

Information in this document is provided solely to enable system and software implementers to use Freescale Semiconductor products. There are no express or implied copyright license granted hereunder to design or fabricate any integrated circuits or integrated circuits based on the information in this document.

Freescale Semiconductor reserves the right to make changes without further notice to any products herein. Freescale Semiconductor makes no warranty, representation or guarantee regarding the suitability of its products for any particular purpose, nor does Freescale Semiconductor assume any liability arising out of the application or use of any product or circuit, and specifically disclaims any and all liability, including without limitation consequential or incidental damages. "Typical" parameters which may be provided in Freescale Semiconductor data sheets and/or specifications can and do vary in different applications and actual performance may vary over time. All operating parameters, including "Typicals" must be validated for each customer application by customer's technical experts. Freescale Semiconductor does not convey any license under its patent rights nor the rights of others. Freescale Semiconductor products are not designed, intended, or authorized for use as components in systems intended for surgical implant into the body, or other applications intended to support or sustain life, or for any other application in which the failure of the Freescale Semiconductor product could create a situation where personal injury or death may occur. Should Buyer purchase or use Freescale Semiconductor products for any such unintended or unauthorized application, Buyer shall indemnify and hold Freescale Semiconductor and its officers, employees, subsidiaries, affiliates, and distributors harmless against all claims, costs, damages, and expenses, and reasonable attorney fees arising out of, directly or indirectly, any claim of personal injury or death associated with such unintended or unauthorized use, even if such claim alleges that Freescale Semiconductor was negligent regarding the design or manufacture of the part.